




1M102 SERIES





KNEES

The 1M102 offers increased safety in stance phase with its alignment and axis geometry.



- Socket rotation adjustment: $\pm 15^\circ$.
- Friction adjustment.
- Extension assist adjustment.
- Optional lock available.



TECHNICAL FEATURES

Item N°.	Max. Patient's Weight	Weight	TH	BH	H	Material	Lock	Flexion	Distal Adapter
1M102	220 lbs/ 100 kg	565 g	184 mm	113 mm	12 mm	Aluminium	No	160 °	Tube Ø 30 mm
1M102V		575 g					Yes		
1M102-P6	275 lbs/ 125 kg	500 g	184 mm	113 mm	12 mm		No		Tube Ø 34 mm
1M102V-P6		525 g					Yes		

ORDER ONE OF THE FOLLOWING TO COMPLETE YOUR ORDER



1K40



1K160-P6



1K160



1K163

ASK A QUESTION

ORDER