

HYTREK

YOUR LIFE. YOUR WAY.

Erin Deegan, AK, HyTrek Knee, RUSH Foot® RAMPAGE LP

A **SINGLE AXIS KNEE** with **HYDRAULIC** regulation of stance phase and swing phase, equipped with a manual lock.

- Strong high-end knee for dynamic transfemoral amputees.
- Navigate everyday life activities on all types of terrain; including slopes and stairs.
- Offers smooth ambulation and the ability to safely descend stairs step over step.

Increased patient dynamics

Optional manual lock

Powerful hydraulic brake

Extension assist

Wide braking range

Maximum patient's weight: 330lbs. (150kg)



STANCE: Adjustment of sensitivity for switching to braking mode.

RESISTANCE: Adjustment of braking power.

LOCK: Manual hydraulic lock switch.

FLEXION: Adjustment of flexion angle.

EXTENSION: Dynamic control of swing phase and terminal impact.



Adjustment wrench (included)



TECHNICAL FEATURES

| | | |
|--------------------------|--------------------------------------|------------------|
| Item Number for Order | 1P131 | 1P131-KD |
| Weight | 2.79lbs. (1265g) | 2.83lbs. (1285g) |
| TH/BH/H (in mm) | 232/177/36 | 222/174/33 |
| Materials | Aluminium, Stainless steel, Titanium | |
| Knee flexion | 120° | |
| Proximal Connector | Pyramid | M36 screw |
| Distal Connector | 34mm Tube | |
| Maximum patient's weight | 330lbs. (150kg) | |
| Activity level | K3-K4 | |
| Water resistance | Weather resistant | |
| Recommended Coding | L5828, L5850, L5925 and L5845 | |

RELATED PRODUCTS

| | | |
|-----------------------|-----------------|----------------------|
| Lamination adapters | 1K179 | 1K179 |
| Adapters | 1K183•1K207-HD | 1K207-HD |
| | 1D41-HD•1K190 | 1D41-HD•1K190 |
| Tubes | 1G01-HD 1G01-HD | |
| | | |
| RUSH Foot® Collection | | 1A600: DynaTrek foot |



1236 West Southern Avenue #101, Tempe, AZ 85282 • 855.450.7300

proteorusa.com #HumanFirst

For additional information on the hyTrek visit: hytrek.proteor.eu

Human First