



| PRODUCT | PRODUCT ID | MODULE | CATEGORY | SIZE | TOE CONFIGURATION |
|----------------|------------|--------|----------|---------|--|
| Agilix | F15 | 00 | 01 - 09 | 22 - 31 | RU: Regular, Universal SR: Sandal, Right* SL: Sandal, Left* <i>*Sandal Toe Available Sizes: 24 - 28: Pacifica & Kinnex 22 - 28: All other models</i> |
| DynAdapt | F10 | | | | |
| Highlander | FS3 | | | | |
| Pacifica LP | FS4 | | | | |
| Sierra | FS1 | | | | |
| Pacifica | FS2 | | | | |
| Highlander MAX | FS3 | H5 | 10 - 12 | 22 - 31 | |
| Kinterra | ROM | 00 | 01 - 07 | 22 - 31 | |
| Kinnex | F14 | N2 | 01 - 07 | 24 - 30 | |

PART NUMBER EXAMPLE: If ordering a Freedom Agilix Foot Module for a 150 lbs. patient with a right-side amputation, moderate activity level, in a size 23 sandal toe, the part number would be: **F15-00-04A23-SR**

| STIFFNESS CATEGORY SELECTION | | | | |
|--|-----------------------|--------------|----------|------|
| PRODUCT | WEIGHT lbs. / kg. | IMPACT LEVEL | | |
| | | LOW | MODERATE | HIGH |
| Agilix* DynAdapt Highlander Pacifica Pacifica LP Sierra | 110 - 115 / 44 - 52 | 1 | 1 | 2 |
| | 116 - 130 / 53 - 59 | 1 | 2 | 3 |
| | 131 - 150 / 60 - 68 | 2 | 3 | 4 |
| | 151 - 170 / 69 - 77 | 3 | 4 | 5 |
| | 171 - 195 / 78 - 88 | 4 | 5 | 6 |
| | 196 - 220 / 89 - 100 | 5 | 6 | 7 |
| | 221 - 255 / 101 - 116 | 6 | 7 | 8 |
| | 256 - 285 / 117 - 130 | 7 | 8 | 9 |
| | 286 - 325 / 131 - 147 | 8 | 9 | |
| Highlander MAX | 286 - 325 / 131 - 147 | 8 | 9 | 10 |
| | 326 - 365 / 148 - 166 | 9 | 10 | 11 |
| | 366 - 420 / 167 - 190 | 10 | | |
| | 421 - 480 / 191 - 218 | 11 | | |
| Kinterra Kinnex | 110 - 115 / 44 - 52 | 1 | 1 | |
| | 116 - 130 / 53 - 59 | 1 | 2 | |
| | 131 - 150 / 60 - 68 | 2 | 3 | |
| | 151 - 170 / 69 - 77 | 3 | 4 | |
| | 171 - 195 / 78 - 88 | 4 | 5 | |
| | 196 - 220 / 89 - 100 | 5 | 6 | |
| | 221 - 255 / 101 - 116 | 6 | 7 | |
| 256 - 275 / 117 - 125 | 7 | | | |

IMPACT LEVELS BY ACTIVITY

| LOW |
|-----------|
| Walking |
| Biking |
| Bowling |
| Camping |
| Dancing |
| Fishing |
| Gardening |
| Golfing |
| Surfing |

| MODERATE |
|---------------|
| Jogging |
| Boxing |
| Hiking |
| Martial Arts |
| Skating |
| Skateboarding |
| Skiing |
| Tennis |
| Volleyball |

| HIGH |
|--------------------|
| Running |
| Competitive Sports |
| Football |
| Power Lifting |
| Snowboarding |
| Soccer |
| Wakeboarding |

**Agilix:*
 - Not recommended for High Impact.
 - Not recommended for 110 - 115lbs. / 44 - 52kg.

FTC -
 _
_
_
_
_
_
_
_
_
_

| PRODUCT | STANDARD FOOT SHELLS | GEN_PROFILE - CAP | SIZE | SKIN TONE | TOE | SIDE |
|--|---|-------------------|----------|---|-------------------|-----------------------------------|
| Agilix DynAdapt Highlander Highlander MAX Pacifica* Pacifica LP Sierra Kinterra | GENERATION 3: Mid Profile, No Cap | 3M - 1 | 22 - 31* | L: Light M: Medium D: Dark | R: Regular | L: Left R: Right |
| | GENERATION 2: Sandal Toe, No Cap | 2M - 1 | 22 - 28* | | S: Sandal | |
| Kinnex | GENERATION 3: Low Profile, No Cap | 3K - 1 | 24 - 30 | | R: Regular | |
| | GENERATION 2: Sandal Toe, No Cap | 2K - 1 | 24 - 28 | | S: Sandal | |

* Pacifica is not available in sizes 22 - 23.

PART NUMBER EXAMPLE: If ordering a Freedom Agilix Foot Shell for a patient with a right-side amputation, in a size 23 sandal toe, and medium skin tone, the part number would be:
FTC-2M-123M4-SR

AVAILABLE UPON REQUEST

| PRODUCT | OTHER FOOT SHELLS | GEN_PROFILE - CAP | SIZE | SKIN TONE | TOE | SIDE |
|--|---|-------------------|----------|---|-------------------|-----------------------------------|
| Agilix DynAdapt Highlander Highlander MAX Pacifica* Pacifica LP Sierra Kinterra | GENERATION 3: Mid Profile, With Cap | 3M - 0 | 22 - 31* | L: Light M: Medium D: Dark | R: Regular | L: Left R: Right |
| | GENERATION 2: High Profile, No Cap | 2H - 1 | | | | |

* Pacifica is not available in sizes 22 - 23.