



# matik



KNEES

## The Matik by Proteor was designed to provide an affordable 4-axis pneumatic knee for patients up to 275 lbs.

Simply and effectively designed to provide high function with low maintenance. The easy swing of the pneumatic unit and the reduced weight of the knee give the patient a heightened sensation of fluidity and comfort. To support safety during swing phase and reduce the risk of a fall, the toe clearance height reaches 25mm, one of the industry's highest for this type of knee. Matik has both flexion and extension adjustments which address terminal impact and excessive heel rise potentials.



- Pneumatic knee with extension spring.
- Secure 4-axis geometry.
- Safety and security through generous toe clearance during swing phase.
- Compact and lightweight (26,21 oz / 743g).
- Adapted for patients up to 275 lbs / 125 kg.

### TECHNICAL FEATURES

|                          |                            |              |
|--------------------------|----------------------------|--------------|
| Item N° for Order        | 1P200                      | 1P200-KD     |
| Weight                   | 743 g                      |              |
| TH/BH/H (In mm)          | 204/175/11,5               | 195/172/8    |
| Materials                | Aluminium, Stainless steel |              |
| Knee flexion             | 170 °                      |              |
| Proximal Connector       | Pyramid                    | Threaded Top |
| Distal Connector         | Pyramid                    |              |
| Maximum patient's weight | 275 lbs / 125 kg           |              |
| Activity level           | 2-3                        |              |

1X998 To avoid spoiling the pyramid during trials



ASK A QUESTION

ORDER