

# RUSH ROVER<sup>®</sup>



## TERRAIN DOMINATION

Introducing the ultra low-profile **RUSH ROVER**. The all-new **RUSH ROVER** features all the same flexibility and durability offered by the RUSH Foot Collection, with a low build height of just over 2"! Same smooth roll-through, same dynamic toe-off, most durable low-profile foot available!

### RUSH ROVER PROFILE INFORMATION

Weight Limit	365lbs/165kg	Heel Rise	3/8" / 10mm	Size	Build Height**
Size	22cm – 30cm	Warranty	36 months	22-24cm	2.125" / 54mm
Foot Shell	Light or Dark	†Suggested Codes	L5981, L5986	25-27cm	2.375" / 60mm
Product Weight*	950g			28-30cm	2.5" / 64mm

† The listing of this code should not be considered a guarantee for coverage or a payment.

\* Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.

\*\*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.



### CATEGORY CHART

Weight	lbs	0 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166	
Low Activity	1	2	3	4	5	6	7	8	9	
Moderate Activity	1	2	3	4	5	6	7	8	9	
High Activity	1	2	3	4	5	6	7	8	9	
High Impact	2	3	4	5	6	7	8	9	-	

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t You in p i #HumanFirst



Human First