



# RISE AND RUN.



**In the gym, at the park, or on the trail... the EasyRun blade by Proteor USA is the perfect foot to get you moving.**

So, you're not an experienced sprinter? Not to worry! Whether you're exploring new activities or looking to take your daily experience to a new level, the **EasyRun** will help you increase strength, mobility, happiness, and amplify your life. Specifically designed for the "every day user", the **EasyRun** is here to help you **#riseandEASYPUN**.

With its all-terrain, all-weather, curved fiberglass blade design, and anti-slip Vibram Sole... you can take your **EasyRun** workout outside with confidence.

**EasyRun** is available in 6 vibrant designs... allowing you to choose one as unique as your personality! Replacement soles are available in both black and lime green to keep your **EasyRun** and on the move.


**THINK STRONG. BE STRONG. FINISH STRONG.**

## Technical FEATURES

<b>Weight</b>	2.40lbs / 1088 g
<b>Build Height</b>	8.66in / 220 mm
<b>Materials</b>	Fiberglass + rubber
<b>Rotatable Pyramid</b>	Pyramid

## PATIENT WEIGHT AND ACTIVITY

<b>Weight</b>	88-108lbs 40-49kg	110-130lbs 50-59kg	132-152lbs 60-69kg	154-174lbs 70-79kg	176-194lbs 80-88kg	198-218lbs 90-99kg	220-242lbs 100-110kg
<b>Moderate Impact</b>		1	2	3	4	5	6
<b>High Impact</b>	1	2	3	4	5	6	



**AVAILABLE  
IN 6 VIBRANT  
DESIGNS**

**#riseandEASYPUN**



1236 West Southern Avenue #101, Tempe, AZ 85282 • **855.450.7300**

proteorusa.com [f](#) [@](#) [v](#) [in](#) **#HumanFirst**

*Human First*



# Human First





ROTATABLE PYRAMID ADAPTER

CURVED FIBERGLASS BLADE

LOWER BLADE WITH  
OPTIMIZED CURVATURE

ANTI-SLIP SOLE



1236 West Southern Avenue #101, Tempe, AZ 85282 • 855.450.7300  
proteorusa.com     #HumanFirst

Human First