

SMOOTH • FLEXIBLE • RELIABLE • VIRTUALLY INDESTRUCTIBLE



RUSH Chopart®

VIRTUALLY INDESTRUCTIBLE

The **RUSH Chopart®** is specifically designed for people with Chopart amputations as well as lisfranc, pirogoff, boyd and partial foot prostheses. A dynamic design, providing increased reliability and eliminating the breakage commonly seen with carbon Chopart plates. The unique glass composite is engineered to be more flexible and durable than any other Chopart available today.



RUSH Chopart PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8" / 10mm
Size	22cm - 30cm	Warranty	24 Months
Foot Shell	Light or Dark	†Suggested Codes	L5976
Product Weight*	380g	Build Height**	1.125" / 29mm

† The listing of this code should not be considered a guarantee for coverage or a payment.

* Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.

**Build height based on size 26cm, Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.

CATEGORY CHART

Weight	0 - 70	71 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
lbs	0 - 32	33 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166
kg										
Low Activity	0	1	2	3	4	5	6	7	8	9
Moderate Activity	0	1	2	3	4	5	6	7	8	9
High Activity	0	1	2	3	4	5	6	7	8	9
High Impact	1	2	3	4	5	6	7	8	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t You Tube in p #HumanFirst



Human First