

DYNAMIC · WATERPROOF · SLIP-PROOF · COMFORTABLE

RUSH FOOT[®] H2O[®] COLLECTION



INCOMPARABLE GRIP. NO FOOT SHELL REQUIRED.

The **RUSH Foot H2O[®] Collection** is the ideal foot for your shower leg or various, non-invasive* water sport activities. Available in the tried-and-true **RUSH Foot[®] HiPro, ROGUE 2, RAMPAGE, and RAMPAGE LP** models and featuring an integrated Vibram[®] sole and custom molded rubber toe wrap, the **RUSH Foot[®] H2O Collection** offers the same great durability and comfort of all the other RUSH Feet, but now adds superior traction on slick surfaces.

From the shower to the pool, you won't find a foot that provides this amount of grip and added confidence on slippery surfaces.

RUSH HiPro H2O PROFILE INFORMATION

Weight Limit	365lbs/165kg	Size	Build Height**
Size	22cm - 29cm	22 - 24cm	6" / 152mm
Foot Shell	No foot shell required	25 - 27cm	6.25" / 159mm
Product Weight†	765g	28 - 29cm	6.875" / 175mm
Warranty	36 Months	†Suggested Codes	L5981

RUSH RAMPAGE LP H2O PROFILE INFORMATION

Weight Limit	365lbs/165kg	Size	Build Height**
Size	22cm - 29cm	22 - 24cm	4.125" / 105mm
Foot Shell	No foot shell required	25 - 27cm	4.125" / 105mm
Product Weight†	650g	28 - 29cm	4.325" / 111mm
Warranty	36 Months	†Suggested Codes	L5981 L5986

RUSH ROGUE 2 H2O PROFILE INFORMATION

Weight Limit	365lbs/165kg	Size	Build Height**
Size	22cm - 29cm	22 - 24cm	5.625" / 142mm
Foot Shell	No foot shell required	25 - 27cm	6" / 152mm
Product Weight†	883g	28 - 29cm	6.5" / 165mm
Warranty	36 Months	***Codes	L5987 L5984



RUSH RAMPAGE H2O PROFILE INFORMATION

Weight Limit	365lbs./165kg	Size	Build Height**
Size	22cm - 29cm	22 - 24cm	6.25" / 159mm
Foot Shell	No foot shell required	25 - 27cm	6.5" / 165mm
Product Weight†	695g	28 - 29cm	7" / 178mm
Warranty	36 Months	†Suggested Codes	L5987 L5986

CATEGORY CHART

Weight	lbs	0 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166
Low Activity	1	2	3	4	5	6	7	8	9	
Moderate Activity	1	2	3	4	5	6	7	8	9	
High Activity	1 (2)	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	
High Impact	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	- (-)	

† The listing of this code should not be considered a guarantee for coverage or a payment.

* Product weight based on size 26cm, Cat 4 foot module.

**Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module.

***PDAC Verified

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability. (RAMPAGE LP)

*The H2O collection feet are for water activities or activities that you typically would do barefoot (such as showering, swimming, surfing, karate, yoga, etc). They are NOT intended for aggressive terrain or extreme element exposure. Excessive non-water related exposure (such as pavement, rocks, concrete, etc) can result in premature composite damage and will void the foot's warranty.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t YouTube in p #HumanFirst

