



MEET THE NEW **KID** ON THE BLOCK

We all know kids can be tough on prosthetic feet; we've got the solution for you! Introducing the ultra low-profile **RUSH Kid**. The all-new **RUSH Kid** features all the same flexibility and durability offered by the RUSH Foot Collection, with a low build height of just 2"! Same smooth roll-through, same dynamic toe-off, most durable pediatric foot available!

RUSH Kid PROFILE INFORMATION

Weight Limit	135lbs / 61kg	Heel Rise	3/8" / 10mm
Size	19cm – 21cm	†Suggested Codes	L5981, L5986
Foot Shell	Light or Dark	Warranty	18 months
Product Weight*	505g	Build Height**	2.0" / 51mm

† The listing of this code should not be considered a guarantee for coverage or a payment.

* Product weight based on size 20cm, Cat 4 foot module with spectra sock and foot shell.

**Build height based on size 20cm, Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.



CATEGORY CHART

Weight	lbs	0 - 30	31 - 45	46 - 60	61 - 75	76 - 90	91 - 105	106 - 120	121 - 135
	kg	0 - 14	15 - 20	21 - 27	28 - 34	35 - 41	42 - 48	49 - 54	55 - 61
Low Activity		1P	2P	3P	4P	5P	6P	7P	8P
Moderate Activity		1P	2P	3P	4P	5P	6P	7P	8P
High Activity		1P	2P	3P	4P	5P	6P	7P	8P
High Impact		2P	3P	4P	5P	6P	7P	8P	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t You in p #HumanFirst



HumanFirst