



Product Code	
Kid	= KID
Rover	= ROV
Chopart	= CHO
HiPro	= HiP
RAMPAGE LP	= RAMLP
RAMPAGE	= RAM
**** RAMPAGE H2O	= H2RAM
**** HiPro H2O	= H2H
**** RAMPAGE LP H2O	= H2RAMLP
**** ROGUE 2 H2O	= H2R2
**** ROGUE 2	= ROG2
* Foot Shell	= FS
** Spectra Sock	= SS
EVAQ8 RAMPAGE	= EVRAM
EVAQ8 ROGUE 2	= EVQR2
EVAQ8 HiPro	= EVQH
EVAQ8 RAMPAGE LP	= EVQRAMLP
EVAQ8	= EVRB
Rebuild Kit	= EVRB
*** EVAQ8 Release Valve	= EVRV

Foot Size	
19 cm*	= 19
20 cm*	= 20
21 cm*	= 21
22 cm	= 22
23 cm	= 23
24 cm	= 24
25 cm	= 25
26 cm	= 26
27 cm	= 27
28 cm	= 28
29 cm	= 29
30 cm**	= 30

* Kid only
** ROVER + Chopart only

Kid CATEGORY CHART

Weight	0 - 30		31 - 45		46 - 60		61 - 75		76 - 90		91 - 105		106 - 120		121 - 135	
	lbs	kg	0 - 14	15 - 20	21 - 27	28 - 34	35 - 41	42 - 48	49 - 54	55 - 61	62 - 68	69 - 75	76 - 82	83 - 89	90 - 96	97 - 103
Low Activity	1P		2P		3P		4P		5P		6P		7P		8P	
Moderate Activity	1P		2P		3P		4P		5P		6P		7P		8P	
High Activity	1P		2P		3P		4P		5P		6P		7P		8P	
High Impact	2P		3P		4P		5P		6P		7P		8P		-	

Foot Shell	
Left Light	= L1
Left Dark	= L3
Right Light	= R1
Right Dark	= R3
H2O	= X

RAMPAGE LP CATEGORY CHART

Weight	0 - 105		106 - 140		141 - 175		176 - 210		211 - 245		246 - 280		281 - 315		316 - 350		351 - 365	
	lbs	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166	167 - 182	183 - 198	199 - 214	215 - 230	231 - 246	247 - 262	263 - 278
Low Activity	1		2		3		4		5		6		7		8		9	
Moderate Activity	1		2		3		4		5		6		7		8		9	
High Activity	2		3		4		5		6		7		8		9		-	
High Impact	3		4		5		6		7		8		9		-		-	

RAMPAGE, ROGUE, HiPro, & ROVER CATEGORY CHART

Weight	0 - 105		106 - 140		141 - 175		176 - 210		211 - 245		246 - 280		281 - 315		316 - 350		351 - 365	
	lbs	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166	167 - 182	183 - 198	199 - 214	215 - 230	231 - 246	247 - 262	263 - 278
Low Activity	1		2		3		4		5		6		7		8		9	
Moderate Activity	1		2		3		4		5		6		7		8		9	
High Activity	1		2		3		4		5		6		7		8		9	
High Impact	2		3		4		5		6		7		8		9		-	

CHOPART CATEGORY CHART

Weight	0 - 70		71 - 105		106 - 140		141 - 175		176 - 210		211 - 245		246 - 280		281 - 315		316 - 350		351 - 365	
	lbs	kg	0 - 32	33 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166	167 - 182	183 - 198	199 - 214	215 - 230	231 - 246	247 - 262	263 - 278	
Low Activity	0		1		2		3		4		5		6		7		8		9	
Moderate Activity	0		1		2		3		4		5		6		7		8		9	
High Activity	0		1		2		3		4		5		6		7		8		9	
High Impact	1		2		3		4		5		6		7		8		9		-	

NOTES

Increasing category numbers indicates increased keel response.
FOR BILATERAL: please increase one category for added stability.

* Customers ordering foot shell only will indicate foot shell number only.

** Customers ordering spectra sock only will indicate spectra sock number only.

*** RUSH EVAQ8 comes with a standard installation kit. Valve not included, but available for additional cost.

**** RUSH H2O does not use foot shell. Please use X in the part number.

PART NUMBER EXAMPLE

If ordering a RUSH Foot HiPro® in size 28cm for a 225 lb user with a moderate activity level, a right-side amputation, and preferring a dark foot shell, the part number would be: **HiP-28-5-R3**

