

KNEES

MPK



L-CODES IN RED = 

	QUATTRO		PLIÉ® 3		ALLUX 2		SYMPHONY				HYTREK		EASYRIDE		MATIK		1M SERIES				1M 102 SERIES				1M 112 SERIES							
Item Number for Order	QNX0010	QNX0011	KS3-00-KNEE1-KT	KD3-00-KNEE2-KT	NE-Z41	NE-Z41SH	NK-6	NK-6+L	NK-6 SH	NK-6 SH+L	1P131	IP131-KD	1P650	1P200	1P200-KD	1M01	1M10	1M01-HD	1M10-HD	1M102	1M102-V	1M102-HD	1M102V-HD	1M112	1M113							
L-Codes	L5856, L5828, L5848, L5845, L5850, L5925		L5828, L5856, L5845, L5848, L5850		K1014, L5856, L5845, L5848		NK-6/SH+L: L5814, L5845, L5848, L5850, L5925 NK-6/SH: L5814, L5845, L5848, L5850				L5828, L5845, L5850, L5925		L5830, L5845, L5848, L5850, L5925, L5930		L5840, L5850		L5811				L5611, L5850, L5925				L5611, L5850							
Product Weight	1655g		1235g		1243g		1510g		1520g		920g	970g	960g	1010g	1265g	1285g	1264g		743g		345g	355g	375g	565g	575g	500g	525g	320g	340g			
Build Height	8.5" / 216mm		9.3" / 235mm		8.8" / 223mm		11.5" / 295mm		11.3" / 287mm		7.8" / 197mm		7.5" / 191mm		7.0" / 177mm		6.9" / 174mm		9.0" / 227mm - 10.0" / 260mm		6.9" / 175mm		6.8" / 172mm		2.3" / 58mm				4.5" / 113mm		3.0" / 75mm	
Materials	Aluminum, SS, Titanium		Aluminum, SS, Titanium		Carbon Fiber, SS		Aluminum, Titanium				Aluminum, SS, Titanium		Aluminum Alloy		Aluminum, SS		Aluminum, Steel				Aluminum				Aluminum							
Knee Flexion	135°		125°		117°		180°		170°		120°		130°		170°		130°				160°				165°							
Proximal Connector	Pyramid		Threaded Top		Pyramid		Threaded Top		Pyramid		Threaded Top		Pyramid		Threaded Top		Pyramid				Threaded Top				Pyramid							
Distal Connector	Pyramid		Pyramid		Pyramid		Pyramid				Pyramid		Pyramid		Pyramid		Pyramid				Pyramid				Pyramid							
Max Patient Weight	300lbs / 136kg		275lbs / 125kg moderate activity 220lbs / 100kg high activity		275lbs / 125kg (K3)		220lbs / 100kg				275lbs / 125kg		330lbs / 150kg		275lbs / 125kg		275lbs / 125kg		220lbs / 100kg		275lbs / 125kg		220lbs / 100kg		275lbs / 215kg		99lbs / 45kg		176lbs / 80kg			
Activity Level	K3 - K4		K3 - K4		K3		K3				K3 - K4		K3 - K4		K3		K1				K1 - K2				K1 - K2							
Water Resistance	IP67		IP67		IP44		-				Weather Resistant		Waterproof		-		-				-				-							

FEET



	RUSH ROGUE® 2	RUSH RAMPAGE®	RUSH RAMPAGE® LP	RUSH HI PRO	RUSH ROVER®	RUSH KID®	RUSH CHOPART	EASYRUN
Item Number for Order	ROG2	RAM	RAMLP	HIP	ROV	KID	CHO	1A920
L-Codes	L5987, L5984	L5987, L5986	L5981, L5986	L5981	L5981	L5981, L5986	L5976	L5980
Product Weight*	1039g	850g	805g	920g	950g	505g	380g	1088g
Max Patient Weight	365lbs / 166kg	365lbs / 166kg	365lbs / 166kg	365lbs / 166kg	365lbs / 166kg	135lbs / 61kg	365lbs / 166kg	88lbs / 40kg - 242lbs / 110kg
Build Height**	22-24cm = 6.1" / 155mm	22-24cm = 6.6" / 168mm	22-24cm = 4.5" / 114mm	22-24cm =				