



EXPECT MORE

As the shortest MPK on the market, PROTEOR QUATTRO[™] Microprocessor Knee allows patients to choose from a wider selection of ankles/feet giving them more of what they need to live the life they want.

PROTEOR QUATTRO AND ANKLE/FOOT COMBINED BUILD HEIGHT

SIZE	22-25cm	26-28cm	29-31cm
KINTERRA [®]	13.15"	13.17"	13.31"
SHOCKWAVE ^{®*}	15.44"	15.94"	16.38"
AGILIX	14.03"	14.75"	15.30"
HIGHLANDER [®]	12.30"	12.84"	13.14"
Pacifica ^{®**}	11.05"	11.10"	11.25"
Pacifica [®] LP	10.40"	10.64"	10.79"
SIERRA [®]	13.81"	14.70"	15.26"
DynAdapt	14.62"	15.05"	15.50"

SIZE	22-24cm	25-27cm	28-30cm
RUSH ROGUE [®] 2 / EVAQ8	14.63"	15.13"	15.25"
RUSH RAMPAGE [®] / EVAQ8	15.13"	15.63"	15.75"
RUSH HiPro [®] / EVAQ8	15.00"	15.38"	15.63"
RUSH RAMPAGE [®] LP / EVAQ8	13.00"	13.25"	13.25"
RUSH Rover [®]	11.63"	10.86"	11.00"

*SHOCKWAVE 23-30cm / **Pacifica 24-31cm

Build Heights Based On:

- PROTEOR QUATTRO Microprocessor Knee: 8.5"
- Ankle/Foot: Size 23cm/26cm/29cm, Cat 4 foot module, foot shell, Spectra[™] sock, and 10mm heel height

Final Build Height Dependent On Socket and Componentry



Experience more with QUATTRO, scan now.