

ENGINEERED TO EXCEED EXPECTATIONS



"I live an active lifestyle and the PROTEOR QUATTRO gives me the edge to do whatever I want to do. I have successfully returned to enjoying all my outdoor activities including riding my bicycle. QUATTRO's cycling activity is one of my favorite features because the knee supports my weight when I step off my bike, instead of collapsing like other knees. The QUATTRO has truly been a game changer for me.



Jim W. | Retired Police Officer

PROTEOR QUATTRO MPK and ShockWave® Foot



The QUATTRO is by far my favorite knee I have used in 24 years as an active amputee. Features such as the intuitive stance lock and obstacle assist make it easy to traverse complex terrain and make me feel safer and more comfortable with my daily life. Plus, the app makes making minor adjustments a breeze and the buttons on the knee allow me to easily lock the knee for tasks such as putting on shoes or pants.

Erin DC. | Senior Accountant

PROTEOR QUATTRO MPK and RUSH® ROGUE 2 Foot



QUATTRO provides me RELIABILITY. I never worry about this knee letting me down. My favorite feature is the ability to make adjustments on the fly through the app. I enjoy the custom golf and cycling activities and the step counter. I also really like that this knee lines up well when I sit, most people are unable to detect that I have a prosthetic knee/leg. This is a knee for all! Anyone with an active to a more sedentary lifestyle can benefit from wearing a QUATTRO, it provides the best mobility and I can always trust and rely on it.



"Dez" Desmond | Elite Athlete

PROTEOR QUATTRO MPK and ShockWave® Foot

TECH SPECS

Knee Flexion: 135°

Build Height: 8.5 in (216 mm)

User Weight Limit: 300 lbs /136 kg

Product Weight: 3.65 lbs (1655 g)

Battery Life: 2-3 days

Warranty: 36 months

Part Numbers:

QUATTRO Knee: QNX0010 (Pyramid) / QNX0011 (Threaded Top)



Trial QUATTRO and experience more

ACCESSORIES

Part Description	Part Number	Included w/QUATTRO	Qty	Available/Sold Separately
QUATTRO Battery Wall Charger & Adapters	QNX0601	YES	1	YES
QUATTRO Instructions for Use (IFU)	IFU-03-014	YES	1	YES
External Battery Booster Pack	ACC0010	NO	1	YES

L-code L5856, L5828, L5848, L5845, L5850, L5925



THE PROTEOR PROMISE

PROTEOR is a leading prosthetic and orthotic device company committed to collaborating with clinicians to spark real-world innovations. We are inspired by the people who use our products. Seeing their everyday achievements motivate us to make more things possible. Our extensive portfolio includes technological breakthroughs that allow the utmost mobility, maneuverability, adaptability, and durability. PROTEOR's process for developing advancements and community support is based on a collaborative approach with clinicians and the people they care for. When we all work together, it's amazing how much patients can achieve.

US Headquarters | 1236 West Southern Ave, Suite 101, Tempe, AZ 85282 855.450.7300 www.us.proteor.com

CONNECT WITH US



© 2024 by PROTEOR | All rights reserved. PROTEOR QUATTRO is a trademark of PROTEOR.



QUATTRO Brochure 1/24

PROTEOR
QUATTRO[™]
MICROPROCESSOR KNEE

EXPECT MORE




ENGINEERED
TO EXCEED
EXPECTATIONS



PROTEOR
QUATTRO[™]
 MICROPROCESSOR KNEE

EXPECT MORE FROM YOUR KNEE


The PROTEOR QUATTRO Microprocessor Knee is designed to be more intuitive, more personalized, and more adaptable. With patented technology, advanced customization options, extended battery life, and simple maintenance, QUATTRO gives you more of exactly what you need to live the life you want.




More intuitive for elevated security and dependability



More personalized for a custom experience



More adaptable for fewer interruptions



EFFORTLESS TRANSITIONS

More intuitive

- QUATTRO's patented Control Technology valve allows infinite levels of resistance for **effortless transitions** between activities.
- Three hyper-responsive sensors read body movements to react to the natural gait of each patient.
- Stumble recovery reacts to each situation by determining the appropriate resistance to **recover more naturally**.



RECOVER MORE NATURALLY



REMOTE FIRMWARE UPDATES



20 ACTIVITIES



More personalized

- With up to **20 activity settings**, allowing you to customize movement for specific activities.
- A **mobile app** provides a simple solution for patients and prosthetists to adjust modes, program, and calibrate settings conveniently from your phone.
- As the shortest MPK on the market, you can choose from a wider selection of ankles and feet.

CLINICIAN AND PATIENT
MOBILE APP



Available on the
App Store

GET IT ON
Google Play



TWO DAY BATTERY LIFE + BATTERY BOOSTER PACK

More adaptable

- A battery life of **two to three days** plus an available Booster Pack with additional days of battery.
- Firmware updates received remotely and **service is not required for three years**.
- **Remote transmission** of 3-months of cadence and usage data.

