

## TECHNICAL SPECIFICATIONS

Weight		77-100 lbs 35-45 kg		101-121 lbs 46-55 kg		122-148 lbs 56-67 kg		149-172 lbs 68-78 kg		173-199 lbs 79-90 kg		200-242 lbs 91-110 kg	
Activity level		1/2	3	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3
Audacieuse HLO1_	> 24 cm*			XS	S	S	M	M	L	L	XL	XL	XL
Intrepide HLO2_	> 19 cm*	PXS	PXS	PXS	PS	PS	PM	PM	PL	PL	PXL	PXL	PXL

Activity level: 1 = Low intensity / 2 = Medium intensity / 3 = Very high intensity (sprint)



**\*Minimum clearance:** distance between the ground and the bottom of the socket on the walking prosthesis with shoe.

**FOR BILATERAL PATIENTS:** Increase by one category for added stability.

# Dare to dream



### BORIS, TRAIL RUNNING

*"This blade has the enormous advantage of being versatile. Whether you're a non-athlete, a regular or an elite athlete, you'll enjoy running with it on all types of terrain."*

### SARAH, HIKING AND ROADS

*"Don't be afraid to challenge yourself, to have no limits. When you try, you realize that you're capable of so much!"*



# HOPPER




# HOPPER

**MORE VERSATILE THAN A TYPICAL RUNNING BLADE.  
MORE DYNAMIC THAN A SPORT FOOT.**



US Headquarters | 1236 West Southern Ave, Suite 101,  
Tempe, AZ 85282 855.450.7300 us.proteor.com



EXCLUSIVELY DISTRIBUTED BY  PROTEOR

## HOPPER BLADE

Engineered with the **highest quality materials** for the utmost **versatility, durability** and **comfort**.

### THE HOPPER BLADE - elevated by AIRBUS™

- **Lightweight, ultra-flexible, dynamic carbon fiber** blade
  - Upcycled from **higher-quality** aerospace-grade material
- **Five levels of stiffness** for diverse types/levels of activity
  - **Returns energy** smoothly and comfortably
- **C-shape** allows **heel-to-toe movement** and **rolling motion**
  - Like a **natural running movement**
- **Built to last**, maintains shape over time



### THE ADAPTOR

**Easy to adjust** with a universal top adaptor

- Both male and female included

### THE OUTSOLE - enhanced by SALOMON™

- Patented sole specifically designed to **grip like a shoe**
  - Optimal for **all terrains**
- Enhances **the step through**, providing a **natural stride** feeling
- Maximizes ground contact surface, **improving stability**
- **Interchangeable**, easy to remove/clean/change

Built for versatile use, on flat or uneven terrain

Used for:

• Running • Dancing • Hiking • Badminton • Tennis • Many other sports



## THE FOUNDER'S WORDS

"My dream was to be able to run with my children, but I was confronted with **different issues**.

This raised the question: How to make sport accessible to amputees?"

Jérôme,  
co-founder of Hopper

## TWO MODELS, DIFFERENT CLEARANCES

### AUDACIEUSE

*Audacious*

Since 2021, the pioneering **Audacieuse** has been helping people around the world rediscover the exhilarating joy of sporting activities.

**Minimum clearance: 24 cm.**

### INTREPID

*Intrepid*

The **Intrepide**, recently developed upon request from users (with feedback from prosthetists) is designed for adventurers of all ages who want to push limits with boundless energy.

**Minimum clearance: 19 cm.**



## OUR MISSION

Sport isn't reserved for elite athletes, nor is it just a casual hobby. It is essential for health, well-being, and social inclusion. **Our mission is to make sport accessible to everyone, regardless of performance or financial considerations.**