


ORDERING GUIDE



CARBON FIBER / ENDURACORE / FUSIONFLEX FEET



Product	ID	–	Category	A/P	Size (cm)*	–	Toe Configuration
Agilix	F15–00	–	01–09	A	22–31	–	RU / SR / SL
Agilix EVAQ8	F15–V2	–	01–09	A	22–31	–	RU / SR / SL
ETHOS LP (pediatrics)	LP3-00	–	01-04	P	19-21	–	RU
ETHOS LP	LP3-00	–	01-09	A	22-30	–	SR / SL
ETHOS LP EVAQ8	LP3-V3	–	01-09	A	26-30	–	SR / SL
DynAdapt	F10–00	–	01–09	A	22–31	–	RU / SR / SL
DynAdapt EVAQ8	F10–V2	–	01–09	A	22–31	–	RU / SR / SL
Highlander	FS3–00	–	01–09	A	22–31	–	RU / SR / SL
Highlander EVAQ8	FS3–V2	–	01–09	A	22–31	–	RU / SR / SL
Highlander MAX	FS3–H5	–	10–12	A	26–31	–	RU / SR / SL
Highlander MAX EVAQ8	FS3–V2	–	10–12	A	26–31	–	RU / SR / SL
Kinnex 2.0	F14–N2	–	01–07	A	24–30	–	RU / SR / SL
Kinterra	RM3–00	–	01–08	A	22–30	–	SR / SL
Kinterra EVAQ8	RM3–V2	–	01–08	A	22–30	–	SR / SL
Pacifica	FS2–00	–	01–09	A	24–31	–	RU / SR / SL
Pacifica LP	FS4–00	–	01–09	A	22–31	–	RU / SR / SL
Sierra	FS1–00	–	01–09	A	22–31	–	RU / SR / SL
Sierra EVAQ8	FS1–V2	–	01–09	A	22–31	–	RU / SR / SL
ShockWave	F20–S3	–	01–08	A	23–30	–	RU / SR / SL
ShockWave EVAQ8	F20–V2	–	01–08	A	23–30	–	RU / SR / SL

*For carbon fiber feet, sandal toe only available for sizes up to 28cm, with the exception of ETHOS LP modules

RU: Regular Universal
SR: Sandal Right
SL: Sandal Left

Ordering and Selection Guide – Hopper

Category selection according to weight and activity													
Weight		77–100 lbs 35–45 kg		101–121 lbs 46–55 kg		122–148 lbs 56–55 kg		149–172 lbs 68–78 kg		173–199 lbs 79–90 kg		200–242 lbs 91–110 kg	
Activity	Minimal Clearance*	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3
Audacieuse HL01	240mm 9.45”	–	–	XS	S	S	M	M	L	L	XL	XL	XL
Intrepide HL02	190mm 7.48”	PXS	PXS	PXS	PS	PS	PM	PM	PL	PL	PXL	PXL	PXL

Activity level: 1 = Fitness walking / 2 = Moderate running / 3 = High intensity running

For bilateral patients: Increase by one category for added stability.



*Minimum clearance: distance between the ground and the bottom of the socket on the walking prosthesis with shoe.



Selection Guide for Foot Category based on Weight and Activity Level

Product	Weight (lbs)	Weight (kg)	Low	Mod	High
Agilix* ETHOS LP DynAdapt Highlander Highlander MAX† Pacifica LP ShockWave** Sierra	100–115	44–52	1	1	2
	116–130	53–59	1	2	3
	131–150	60–68	2	3	4
	151–170	69–77	3	4	5
	171–195	78–88	4	5	6
	196–220	89–100	5	6	7
	221–255	101–116	6	7	8
	256–285	117–130	7	8	9
	286–325	131–147	8	9	10
	326–365	148–166	9	10	11
Kinnex 2.0	366–420	167–191	10	–	–
	421–480	192–218	11	–	–
	481–500	219–226	12	–	–
	100–115	44–52	1	1	–
	116–130	53–59	1	2	–
	131–150	60–68	2	3	–
	151–170	69–77	3	4	–
	171–195	78–88	4	5	–
	196–220	89–100	5	6	–
	221–255	101–116	6	7	–
Kinterra	256–275	117–130	7	–	–
	100–115	44–52	1	1	–
	116–130	53–59	1	2	–
	131–150	60–68	2	3	–
	151–170	69–77	3	4	–
	171–195	78–88	4	5	–
	196–220	89–100	5	6	–
	221–255	101–116	6	7	–
	256–290	117–132	7	8	–
	291–330	133–150	8	–	–
PEDS					
ETHOS LP	58	>26	–	–	1P
	9-70	27-36	–	–	2P
	59-79	37-45	–	–	3P
	101-132	46-60	–	–	4P

*Agilix is recommended for only low to moderate impact level.

**ShockWave only goes up to a category 8

†Category 10–12 only available with Highlander MAX

NOTE: Increasing category numbers indicate increased keel response.

Ordering Guide – Gery and DynaStep



Product	Item Number	–	Side	Size (cm)	Color
Gery*	1A200	–	D / G	22–29	H / BH
DynaStep	1A101	–	D / G	22–29	–
DynaStep Foot + Multi-axis Ankle Kit	1E201	–	D / G	22–29	–


*Gery is not recommended for bilateral patients

D = Right
G = Left

H = Light
BH = Dark



FIBERGLASS FEET



Product	ID	–	Size (cm)*	–	Category	–	Toe Configuration
RUSH Chopart	CHO	–	19–30	–	1–9	–	RU
RUSH HiPro	HIP	–	22–29	–	1–9	–	RU, SR, SL
RUSH HiPro EVAQ8	EVQH	–	22–29	–	1–9	–	RU, SR, SL
RUSH HiPro H2O	H2H	–	22–29	–	1–9	–	RU
RUSH RAMPAGE	RAM	–	22–29	–	1–9	–	RU, SR, SL
RUSH RAMPAGE EVAQ8	EVRAM	–	22–29	–	1–9	–	RU, SR, SL
RUSH RAMPAGE H2O	H2RAM	–	22–29	–	1–9	–	RU
RUSH RAMPAGE LP	RAMLP	–	22–29	–	1–9	–	RU, SR, SL
RUSH RAMPAGE LP EVAQ8	EVQRAML	–	22–29	–	1–9	–	RU, SR, SL
RUSH RAMPAGE LP H2O	H2RAMLP	–	22–29	–	1–9	–	RU
RUSH ROGUE 2	ROG2	–	22–29	–	1–9	–	RU, SR, SL
RUSH ROGUE 2 EVAQ8	EVQR2	–	22–29	–	1–9	–	RU, SR, SL
RUSH ROGUE 2 H2O	H2R2	–	22–29	–	1–9	–	RU

RU: Regular Universal
SR: Sandal Right
SL: Sandal Left



SCAN FOR FOOT SHELL INFORMATION

Selection Guide for Foot Category based on Weight and Activity Level

Product	Weight (lbs)	Weight (kg)	Low	Mod	High
RUSH Chopart* RUSH HiPro RUSH RAMPAGE RUSH ROGUE 2	0–105	0–48	1	1	2
	106–140	49–64	2	2	3
	141–175	65–79	3	3	4
	176–210	80–95	4	4	5
	211–245	96–111	5	5	6
	246–280	112–127	6	6	7
	281–315	128–143	7	7	8
	316–350	144–159	8	8	9
	351–365	160–166	9	9	–
RUSH RAMPAGE LP	0–105	0–48	1	2	3
	106–140	49–64	2	3	4
	141–175	65–79	3	4	5
	176–210	80–95	4	5	6
	211–245	96–111	5	6	7
	246–280	112–127	6	7	8
	281–315	128–143	7	8	9
	316–350	144–159	8	9	9
	351–365	160–166	9	–	–
PEDS					
RUSH Chopart	0–30	0–14	1P	1P	2P
	31–45	15–20	2P	2P	3P
	46–60	21–27	3P	3P	4P
	61–75	28–34	4P	4P	5P
	76–90	35–41	5P	5P	6P
	91–105	42–48	6P	6P	7P
	106–120	49–54	7P	7P	8P
	121–135	55–61	8P	8P	9P

*For RUSH Chopart foot sizes 22cm and above only. Reference the PEDS Selection Chart for sizes below 22cm.

NOTE: Increasing category numbers indicate increased keel response.

FOR BILATERAL PATIENTS: Increase by one category for added stability.

CONNECT WITH US



Global Headquarters | 1236 West Southern Ave, Suite 101, Tempe, AZ 85282 | Phone: 855.450.7300 [us.proteor.com](https://www.us.proteor.com)

© 2025 by PROTEOR | All rights reserved.

FOOT SHELLS & SPECTRA SOCKS

FEET

PROTEOR offers a variety of cosmetic foot shells and Spectra Socks.

Please use the following guide to determine which foot shells and socks are available for each foot product. Not all Foot Modules are available in Sandal Toe; ensure the corresponding Sandal Toe foot module is also ordered (CAPs are not available for sandal toe products).

Carbon Fiber / EnduraCore / FusionFlex Modules		___-___-__	___	___	___-__	___
Product	Foot Shell Type	Product ID	Size	Shell Color L: Light / M: Medium D: Dark / B: Jet Black	Toe Configuration	Side L: Left R: Right
PROTEOR Agilix* DynAdapt* Highlander* Highlander MAX* Pacifica LP Sierra*	Regular Toe, with Cap	FTC-3M-0	22-31	L / M / D	4-R	L / R
	Regular Toe, No Cap	FTC-3M-1	22-31	L / M / D / B	4-R	L / R
	Sandal Toe, No Cap	FTC-2M-1	22-28	L / M / D / B	4-S	L / R
ETHOS LP PEDS	Regular Toe, No Cap	FTC-2M-1	19-21	L / M / D / B	4-R	L / R
ETHOS LP*	Sandal Toe, No Cap	FTC-2K-1	22-30	L / M / D / B	4-S	L / R
ShockWave*	Regular Toe, No Cap	FTC-3F-1	23-30	L / M / D / B	4-R	L / R
	Sandal Toe, No Cap	FTC-2F-1	23-28	L / M / D / B	4-S	L / R
PROTEOR Kinnex 2.0	Regular Toe, No Cap	FTC-3K-1	24-30	L / M / D / B	4-R	L / R
	Sandal Toe, No Cap	FTC-2K-1	24-28	L / M / D / B	4-S	L / R
Kinterra*	Sandal Toe, No Cap	FTC-2K-1	22-30	L / M / D / B	4-S	L / R

*Includes EVAQ8 foot modules

**Foot shell compatibility for the older Kinterra 2.0 module: FTC-3M and FTC-2M

FOOT SHELL PART NUMBER EXAMPLE: To order an Agilix Foot Shell in a size 26cm right-side sandal toe, medium skin tone, the part number would be: FTC-2M-126M4-SR



CONTACT US

ORDER

f i y t #HUMANFIRST

© 2025 by PROTEOR | All rights reserved.

us.proteor.com



FOOT SHELLS & SPECTRA SOCKS

FEET

PROTEOR offers a variety of cosmetic foot shells and Spectra Socks.

Please use the following guide to determine which foot shells and socks are available for each foot product. Not all Foot Modules are available in Sandal Toe; ensure the corresponding Sandal Toe foot module is also ordered (CAPs are not available for sandal toe products).

RUSH — Sandal Toe Foot Shell		— — — — —	— —	—	— — —	—
Product	Foot Shell Type	Product ID	Size	Shell Color L: Light / M: Medium D: Dark / B: Jet Black	Toe Configuration	Side L: Left R: Right
RUSH ROGUE 2* RUSH HiPro* RUSH RAMPAGE* RUSH RAMPAGE LP*	Sandal Toe, No Cap	FTC-4R-1	22-29	L / M / D / B	4-S	L / R

RUSH — Regular Toe Foot Shell		— — — — —	— — — — —	— — — — —	—
Product	Foot Shell Type	Product ID	Size	Side L: Left R: Right	Shell Color L: Light D: Dark
RUSH ROGUE 2* RUSH HiPro* RUSH RAMPAGE* RUSH RAMPAGE LP*	Regular Toe, with Cap	FS	22-29	L / R	L / D
RUSH Chopart	Regular Toe, with Cap	FS	22-30	L / R	L / D
RUSH Chopart PEDS	Regular Toe, with Cap	FSK	19-21	L / R	L / D

Spectra Socks



Description	Item N°.
Black Spectra Sock (1 pk)	
22-25cm foot module = 30cm/12"	SO-NPS-200XX-00 (example: size 26cm foot module = SO-NPS-20030-00)
26-28cm foot module = 35cm/14"	
29-30cm foot module = 40cm/16"	
Kinnex White Custom Spectra Sock (1 pk)	
Sizes: 24 - 30cm	SO-F14-000XX-00 (example: size 26cm foot module = SO-F14-00026-00)

*Includes EVAQ8 foot modules

FOOT SHELL PART NUMBER EXAMPLE: To order a RUSH ROGUE 2 Foot Shell in a size 26cm right-side sandal toe, medium skin tone, the part number would be: FTC-4R-126M4-SR



CONTACT US

ORDER

f i y t #HUMANFIRST

© 2025 by PROTEOR | All rights reserved.

us.proteor.com

