



**REQUIRED POWER ADAPTER:
18W+ [9V or 15V]**



USB A



USB MICRO B



USB C



18W+ [9v or 15v]

FR	⚠ Adaptateur secteur requis: 18W+ [9V ou 15V]	FI	⚠ Vaadittu virtälähdeadapteri: 18W+ [9V tai 15V]
DE	⚠ Erforderlicher Netzadapter: 18W+ [9V oder 15V]	CS	⚠ Požadovaný napájecí adaptér: 18W+ [9V nebo 15V]
IT	⚠ Adattatore di alimentazione richiesto: 18W+ [9V o 15V]	SK	⚠ Požadovaný napájací adaptér: 18W+ [9V alebo 15V]
ES	⚠ Adaptador de corriente requerido: 18W+ [9V o 15V]	HR	⚠ Potreban adapter za napajanje: 18W+ [9V ili 15V]
PT	⚠ Adaptador de energia necessário: 18W+ [9V ou 15V]	RU	⚠ Требуемый адаптер питания: 18W+ [9V или 15V]
NL	⚠ Vereiste voedingsadapter: 18W+ [9V of 15V]	UK	⚠ Необхідний адаптер живлення: 18W+ [9V або 15V]
DA	⚠ Krævet strømforsyningsadapter: 18W+ [9V eller 15V]	JA	⚠ 必要な電源アダプター: 18W以上 [9Vまたは15V]
NO	⚠ Påkrevd strømadapter: 18W+ [9V eller 15V]	ZH	⚠ 所需电源适配器: 18W以上 [9V或15V]
SV	⚠ Krävd strömadapter: 18W+ [9V eller 15V]	AR	⚠ بولطم الي ئاياره كذا ل و ح م ل ت ل و ف 15 و أ ت ل و ف 9 ر ك أ ف ط و 18

US: 1236 West Southern Ave, Suite 101 | Tempe, AZ 85282
855.450.7300 | us.proteorusa.com

FRANCE: 6 rue de la Redoute | 21850 | Saint-Apollinaire | France
+33(0)3 80 78 42 42 | fr.proteorusa.com

