



PROTEOR

QUATTRO 2TM

MICROPROCESSOR KNEE

Experience the Evolution.

#HumanFirst



PROTEOR

→ What's New with QUATTRO 2

Not just an upgrade. Next-level features, designed for you.

QUATTRO 2 doesn't just sense your movement... it adapts in real time. Everything that made the original legendary—now stronger, smarter, and more capable, with movement possibilities unique to QUATTRO 2.

A reinforced frame. A more powerful processor.



Stair Ascent Functionality

Natural, intuitive stair climbing powered by a faster processor and reinforced frame.



IP-68 Water Rating

Rain, puddles, the shower—live your life without thinking twice.



150kg / 330lb Weight Limit

A stronger frame expands access to more patients than ever.



10× Faster Microprocessor

A faster, more powerful processor means the QUATTRO 2 responds to your movement in real time, so you can focus on living, not thinking about your knee.



Stronger Frame

Engineered for higher weight capacity and long-term durability.



USB-C Power OPTIONAL

2–3 day battery, one cable, no fuss

And much more.

→ Built Around You

3 Core Promises.

At its core, QUATTRO 2 is built around three commitments.
Motion. Adaptability. Proof.



MOTION

Designed for Natural Movement

A patented single bidirectional valve delivers infinite resistance settings with ultra-responsive speed. Three high-precision sensors drive real-time adaptive adjustments. **The result: smooth, natural, optimized movement.**

ADAPTABILITY

Seamless Personalization

At just 21.6 cm (8.5 in), QUATTRO 2 is the shortest MPK on the market, expanding compatibility with the widest range of feet, ankles, and running blades. With up to 20 pre-set activities, independent stair ascent, stair descent, and ramp resistance, it adapts to whatever your day demands. And with remote firmware updates, your knee stays current—always.

PROOF

Proven by Data, Not Just Claims

In a multicenter, randomized cross-over trial, QUATTRO significantly improved daily life activities compared to the typically worn MPK. Subjects also reported fewer physical limitations with QUATTRO.

→ Clinical Evidence

The data backs every claim.

QUATTRO was compared head-to-head against each subject's currently worn MPK* across real-world daily activities.

We wanted to find out where QUATTRO stood, so we conducted a multicenter, randomized cross-over trial. Our primary objective was to evaluate subjects' ability to achieve personal goals, including level of difficulty. Then we wanted to know their satisfaction and quality of life. Here's what we learned.

4/4

Uneven terrain w/ obstacles

3/3

Walking downhill

2/2

Walking backwards

4/6

Water activities


4/4, for example, means all four users who tested uneven terrain saw.

Key Findings

- **Improved results** in specific daily life activities compared to the typically worn MPK
- **Significantly fewer** physical limitations with QUATTRO than with their typically worn MPK

*RCT, n=17, France, 2023–2024. Protocol approved by national ethics committee (Protocol ID-RCB 2022-A01266-37).

PSFS Scores for Identified Daily Life Activities



*RCT, n=17, France, 2023–2024.

†Number of users that increased their score with QUATTRO compared to typically worn MPK

→ Tech Specs

Just the facts.

All you need to know—now stronger, smarter, and more adaptable than ever.

Core Specifications

Knee Flexion 135°

Build Height 21.6 cm (8.5 in)

Weight Limit 150 kg / 330 lbs **NEW**

Product Weight 1655 g (3.65 lbs)

Battery Life 2–3 Days

Water Rating IP-68 **NEW**

Warranty 36 Months

Optional 2- and 3-year warranty extensions available

Part Numbers


Pyramid QNX0010

Threaded Top QNX0011

Booster Cable (USB-C) ACC0010



L-CODES · PDAC VERIFIED · L5856 · L5828 · L5848 · L5845 · L5850 · L5925





Human First. **Living Proof.**

Shaped by real people and real experiences. QUATTRO 2 was designed with you in mind—built to fit your life and your needs. A connected app puts you at the center of it all, with effortless integration and intuitive controls. When you put people first, innovation follows.


TRIAL NOW

Try QUATTRO Risk-Free.



ORDER NOW

shop.proteorusa.com



info@proteor.com | shop.proteorusa.com | 855-450-7300

[f](#) [YouTube](#) [in](#) [Instagram](#) #FluentInMovement #HumanFirst